TECHNICAL DATA SHEET ⁰

Product Type: | **Bodyjets** Product Code: **E044261**

Product Name: | MIST Bodyjet


Description: | Bodyjet con getto nebulizzato

Bodyjet with nebulized jet


Dimensions: 80x80 mm

Material: ABS


Note: connessione | completamente ispezionabile
Notes: connections G1/2" | completely inspectionable



DIMENSIONS



FLOWRATE



PRODUCT CARE

To avoid damages to the purchased equipment, please observe the following recommendations for use:

- Clean the plumbing before connecting the device.
- Periodically clean the water blade to avoid any limestone deposits that could compromise the good functioning.

Advice: The installation of a filter and water-softener system improves and extends the life of the devices and of the plant itself.

CLEANING OF METAL PARTS:

For an easy cleaning of the metal parts of the product, simply use a soft cloth moistened in water and soap or an anti-limestone detergent.

To avoid spots caused by calcium in water, we advise drying with a cloth after each use.

Advice: Please use only natural soap-based detergents avoiding those with alcohol, hydrochloric acid or phosphoric acid, bleach, chlorine etc. and those containing abrasive powders.

PACKING



6 pcs